


SNAP-LOC SIP-N-PUFF

PARTS DESCRIPTION



REPLACEMENT PARTS


ITEM	PART#	Figure-8 Clamp	32310
Air Tube Straw/10Pack	32698	Figure-8 Clamp, Black	32310B
Air Tube Straw/100Pk	32662	Large Fig.-8, Black	32312
Saliva Trap/5Pack	32661	Multimount Adaptor	32826
1/8" I.D. Clear Vinyl Tube, 5 yards			30088



9450 W Laraway Rd • Frankfort, IL 60423
 Customer 800-843-7234 • Fax 888-479-1515
 www.Therafin.com • info@Therafin.com

MOUNTING SYSTEM TO CHAIR FRAME

1. Use allen wrench provided to loosen **Figure-8 Clamp's** socket head cap screws. Place **Figure-8 Clamp** with Sip-N-Puff system onto chair frame*, and tighten socket head cap screws. The **Figure-8 Clamp** can be turned and adjusted vertically during final adjustment.



***NOTE:** When connecting to a Multimount Chair Frame with a hole pattern of 3/4" or 5/8" on center, use Therafin's **Multimount Adaptor (A)**. LaBac's **Therafin Chin Adaptor (B)** (LaBac part # 11-3J) has a 5/8" on center hole pattern, and is used exclusively with LaBac systems.

2. For final adjustment, bend flexible **Gooseneck** and turn or vertically slide **Figure-8 Clamp** as needed to allow the **Air Tube Straw** to be positioned near the user's mouth. Re-tighten socket head cap screws. Set screws may be tightened, if desired.

MAINTENANCE

To change **Air Tube Straw**, loosen and remove compression fitting nut. Push straw through nut, and cut off compressed end of straw. Pull straw through nut and discard. To insert new **Air Tube Straw** into nut, place bottom end of nut onto hard, flat surface (table top). Push straw end down into hole, through compression flange, until straw makes contact with table top. Retighten nut onto compression fitting.

CLEANING

For external cleaning, wipe with moist cloth. Use mild detergent if desired. For internal cleaning, disconnect from ECU, remove saliva trap, and flush with warm running water. If available, rinse with a disinfecting oral rinse (Listerine, Scope, etc.). Cleaning at least twice a week is recommended, more as needed.